

Dialysis Material	Brand name	Manufacturer	MWCO (kDa)	References
Cellulose esters, mixed	Biotech CE	Spectrum Labs	0.5 to 1000	Biotech CE Data
Cellulose, Regenerated	Dialysis Tubing	Sigma-Aldrich	12.4	Dialysis Tubing Data
Cellulose, Regenerated	Dialysis Tubing	Sigma-Aldrich	14-Dec	Dialysis Tubing Data
Cellulose, Regenerated	Dialysis Tubing, Benzoylated	Sigma-Aldrich	1.2 - 2	Dialysis Tubing, Benzoylated Data
Cellulose, Regenerated	Biotech RC	Spectrum Labs	3.5 - 5, 8 - 10, 20	Biotech RC Data
Cellulose, Regenerated	Servapor HMF	Serva	1 - 25	Servapor HMF Data
Cellulose, Regenerated	Visking	Medicell Membranes	3.5, 7, 14	Visking Data
Cellulose, Regenerated	Cuprophane	Medicell Membranes	5, 10	Cuprophane Data
Cellulose, Regenerated	Ultracel	EMD Millipore	1 to 100	Ultracel Data
Polyamide	Dialysis Membrane	Metrohm	1000	Dialysis Membrane Data
Polyethersulfone PES	Biomax	EMD Millipore	5 to 500	Biomax Data
Polyethersulfone PES	Omega	Pall	1 to 300	Omega Data